


INTRODUCTION

Attachment of polyethylene glycol (PEG) chains to therapeutic biologic agents, a process referred to as PEGylation, prolongs the circulating half-life of the modified protein by slowing proteolytic degradation and by masking it from the immune system. The pharmacodynamics of PEGylated proteins are often evaluated using specific assays for the protein itself. That approach often requires the time consuming and expensive development of an ELISA for the protein of interest. The PEG backbone ELISA allows measurement of the PEG portion of the PEGylated protein and is therefore suitable for assessment of the pharmacodynamics of a range of PEGylated proteins. As shown in the figure below, sensitivity of the assay varies with chain-length and the extent of PEGylation. To determine PEG concentration in unknown samples, the end-user must establish a standard curve using the PEG reagent under investigation.



PRINCIPLE OF THE ASSAY

The assay uses two different monoclonal antibodies that recognize the PEG backbone. One is coated on wells of the microtiter plate and is used for capture; the other is conjugated to horseradish peroxidase (HRP) and is used for detection. Serum samples are diluted at least 10-fold in the provided dilution buffer and incubated alongside reference standards¹ in the microtiter wells for 45-minutes. The wells are subsequently washed. HRP conjugate is added and incubated for 45 minutes. If present, PEG molecules are sandwiched between the capture and detection antibodies. The wells are then washed to remove unbound HRP-conjugate. TMB is added and incubated for 20 minutes. This results in the development of a blue color. Color development is stopped by the addition of Stop Solution, changing the color to yellow. Absorbance is measured at 450 nm. The concentrations of PEG in the samples are derived from a standard curve.

MATERIALS

Materials provided with the kit:

- Anti-PEG coated plate (12 x 8-wells)
- Anti-PEG HRP Stock (lyophilized)
- Test Stock (lyophilized)
- 20x HRP PEG Wash: PEGW50-20, 50 ml
- PEG Diluent: PEGBC50-1, 50 ml
- TMB: TMB11-1, 11 ml
- Stop Solution: SS11-1, 11 ml

Store at -20°C

Store at -20°C

Store at -20°C

Materials required but not provided:

- Pipettors and tips
- Distilled or deionized water
- Polypropylene or glass tubes
- Vortex mixer
- Absorbent paper or paper towels
- Plate incubator/shaker
- Plate washer
- Plate reader capable of measuring at 450 nm.
- Curve fitting software

STORAGE

The test stock, HRP conjugate, and the antibody coated plate should be stored at -20°C. All remaining kit components should be stored at 4°C. The microtiter plate should be kept in a sealed bag with desiccant. Kits will remain stable for six months from the date of purchase provided that the components are stored as described.

GENERAL INSTRUCTIONS/LIMITATIONS

1. Please fully read and thoroughly understand the instructions before using the kit.
2. All reagents should be allowed to reach room temperature (25°C) before use.
3. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
4. Use only the wash solution and dilution buffers provided with the kit. They are specially formulated for measurement of PEG.
5. Kits are validated using plate shakers set at 150 rpm and 25°C. Performance of the assay at lower temperatures and/or mixing speeds will result in lower absorbance values.
6. Optimal results are achieved if at each step, reagents are pipetted into the wells of the microtiter plate within 5 minutes.

¹ The end-user must establish a standard curve using the PEG, or PEGylated reagent, under investigation. Sensitivity of the kit varies with PEG chain length and the extent of PEGylation. The "Test Stock" provided with the kit should only be used to confirm that the assay is working as intended.

WASH SOLUTION

The wash solution is provided as a 20x stock. Prior to use, dilute the contents of the bottle (50 ml) with 950 ml of distilled or deionized water. DO NOT substitute other wash solutions.

DILUENT

The diluent (PEGBC50-1) is formulated for measurement PEG. It is supplied ready to use. DO NOT substitute other buffers.

TEST STANDARDS (OPTIONAL)

The Test Stock provided with the kit consists of BSA to which a single chain of 20 kDa mPEG is attached. It is provided for optional use to confirm that the assay is working as intended. As described in the next section, the end user must establish a standard curve using the PEG reagent under investigation.

1. Reconstitute the Test Stock as described on the vial label.
2. Label 8 microcentrifuge tubes² as 50, 16.67, 5.56, 1.852, 0.617, 0.206, 0.069, and 0.023 ng/ml.
3. In the tube labeled 50 ng/ml prepare the 50 ng/ml standard as detailed on the stock vial label.
4. Dispense 250 μ l of diluent into the remaining tubes.
5. Prepare a 16.67 ng/ml standard by diluting and mixing 125 μ l of the 50 ng/ml standard with 250 μ l of diluent in the tube labeled 16.67 ng/ml.
6. Similarly prepare the remaining standards by three-fold serial dilution.

USER-DEFINED STANDARD CURVE

Because the sensitivity of the ELISA varies with PEG chain length and the extent of PEGylation, the end-user must establish a standard curve using the PEG reagent under investigation. We suggest that an initial dilution series be prepared as follows

1. Prepare a 1000 ng/ml stock of the PEG reagent in diluent, PEGBC50-1.
2. Label 7 microcentrifuge tubes as 200, 40, 8, 1.6, 0.32, 0.64 and 0 ng/ml and dispense 200 μ l of diluent into each.
3. In the tube labeled 200 ng/ml mix 50 μ l of the 1000 ng/ml with 200 μ l of diluent.
4. Similarly prepare the 40 to 0.64 ng/ml standards by five-fold serial dilution.

Test the standards in singlets using one 8-well strip. Use the results as a guide to optimize the standard range. When the high standard is defined (ideally giving an A_{450} of 3 to 3.5), we recommend using three-fold dilutions to prepare lower standards in subsequent experiments.

SAMPLES

In studies at Life Diagnostics, Inc., we performed spike-recovery experiments using rabbit serum. When serum was spiked with 250 or 50 ng/ml of 20 kDa mono mPEG-BSA, respective recoveries of 244 and 39 ng/ml were obtained. Dilutional linearity was observed. To avoid matrix effects serum must be diluted at least ten-fold in diluent PEGBC50-1. The end user must determine optimal dilutions.

HRP CONJUGATE

Approximately 15 minutes before needed, reconstitute the lyophilized HRP conjugate as directed on the vial label and mix gently. Then, dilute as described on the vial label to give the working conjugate solution. The reconstituted conjugate stock should be stored at or below -20°C in a sealed vial if future use is intended.

PROCEDURE

1. Secure the desired number of coated wells in the holder.
2. Dispense 100 μ l of standards and diluted samples into the wells (we recommend testing in duplicate).
3. Incubate on a plate shaker at 150 rpm/25°C for 45-minutes.
4. Aspirate the contents of the microtiter wells and wash the wells five times with 1x wash solution using a plate washer (400 μ l/well).
5. Strike the wells sharply onto absorbent paper, if necessary, to remove all residual wash solution.
6. Add 100 μ l of diluted HRP conjugate into each well.
7. Incubate on a plate shaker at 150 rpm/25°C for 45-minutes.
8. Wash as detailed above.
9. Dispense 100 μ l of TMB into each well.
10. Incubate on a plate shaker at 150 rpm/25°C for 20 minutes.
11. Stop the reaction by adding 100 μ l of Stop Solution to each well.
12. Gently mix. It is important to make sure that all the blue color changes to yellow.
13. Read the optical density at 450 nm with a microtiter plate reader within five minutes³.

² If possible, we recommend that standards and samples be prepared in a blank 96-well plate. This allows for rapid transfer to the assay plate using multi-well pipettors.

³ On certain plate readers the A_{450} value of the high standard may be out of range. If that occurs, absorbance values for all wells may be read at 405 nm instead. Absorbance values will be lower, but this does not affect results.


RESULTS

1. Using curve fitting software, construct a standard curve by plotting absorbance values of the standards versus \log_{10} of the concentration.
2. Fit the standard curve to a four-parameter logistic regression (4PL) equation (x axis = \log_{10} concentration) and determine the concentration of the samples from the standard curve (remember to derive the antilog).
3. Multiply the derived concentration by the dilution factor to determine the actual concentration in the samples.
4. If the A_{450} values of samples fall outside the standard curve, samples should be diluted appropriately and re-tested.

TEST STANDARD CURVE

A standard curve for the test reagent (20 kDa mono-PEGylated BSA) with absorbance at 450nm on the Y-axis and \log_{10} PEG-BSA concentration on the X-axis is shown below.

PEG-BSA (ng/ml)	A_{450}
50	3.245
16.67	2.680
5.556	1.553
1.852	0.864
0.617	0.475
0.206	0.297
0.069	0.226
0.023	0.191



Rev 021225