

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Anti-PEG Agarose***

Catalog No.	9B5-AG
Description	Anti-PEG monoclonal antibody (9B5-6-25-7) coupled to agarose, ~0.6 mg/ml
Applications	Affinity purification of PEGylated proteins
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, 0.09% NaN <sub>3</sub> , pH 7.2
Storage	-2-8°C. DO NOT FREEZE

9B5 Anti-PEG agarose is prepared by irreversibly coupling monoclonal antibody 9B5-6-25-7 to agarose via its FC region. This ensures that the PEG binding sites of the coupled antibody retain functionality. Approximately 0.6 mg of antibody is coupled per ml of settled agarose.

P.O. Box 5205 • West Chester, PA 19380 • USA