

# *Life Diagnostics, Inc.*

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## Acridan Conjugation Service

Catalog No.	SP-AC
Description	Custom acridan conjugation of antibodies for use in SPARCL™ assays.
Antibody Requirement	Antibody must be provided to LDI at a concentration of 1 mg/ml or greater.
Minimum Quantity	1 mg

**Comments.** Acridan conjugated antibodies will be returned to the customer in a proprietary stabilization buffer that we use in our SPARCL™ kits. The conjugate will be shipped on dry ice and should be stored at -80°C. We find that most acridan conjugated antibodies are stable for several years at -80°C and that they withstand several freeze-thaw cycles. For optimum stability, however, we recommend that conjugates be stored as multiple aliquots after the initial thaw. Thaw the conjugate at room temperature.

Rev 030317

P.O. Box 5205 • West Chester PA 19380 • USA