

# *Life Diagnostics, Inc.*

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## SPARCL™ Plate

(12 x 8-well strips)

Catalog No.	SP-PL
Analysis	SPARCL™ assay. CV's ≤ 10%
Coating	Proprietary background reducer.
Storage	-20°C

Comments: For use in SPARCL™ assays. Provided as 12 x 8-well removable strips coated with a proprietary background reducer. Total assay volume should not exceed 100 µl.

Rev 030317

P.O. Box 5205 • West Chester PA 19380 • USA