

# *Life Diagnostics, Inc.*

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## SPARCL™ Conjugate Buffer (25 ml)

Catalog No.	SPB25-1
Description	SPARCL™ Conjugate Buffer
Composition	A proprietary buffer for storage of acridan and HRP conjugate stocks used in SPARCL™ assays.
Storage	2-8°C

**Acridan conjugates.** Conjugates should be prepared as described in the Lumigen SPARCL™ kit instructions. After preparation, dilute the conjugate with 3 volumes of SPB25-1 and stored as aliquots at -80°C (storage at -20°C for 1 week is usually ok but depends on the antibody).

**HRP conjugates.** HRP conjugates, in azide-free buffer, should be diluted with 3 volumes of SPB25-1 and stored as aliquots at -80°C (storage for several weeks at -20°C is usually ok but depends on the antibody).

Rev 030317

P.O. Box 5205 • West Chester PA 19380 • USA