




PEG PROTEINS




MULTI-PEGYLATED BSA, PBSA-00-550, 1 MG

Prepared by reaction of BSA with a 90-fold molar excess of 550 Da mPEG...

SKU: PBSA-00-550

Price: \$500.00

Categories: [PEG Proteins](#), [PEG Reagents](#)



MONO-PEGYLATED BSA (MPEG 20 KDA) BIOTIN CONJUGATE

Prepared by conjugation of mPEG-maleimide (20 kDa) to the single reactive...

SKU: PBSA-01-1

Price: \$600.00 – \$5,500.00

Categories: [PEG Proteins](#), [PEG Reagents](#)



MONO-PEGYLATED BSA (PENTANOIC ACID-PEG, 20 KDA), PBSA-01P, 0.1 MG

Prepared by conjugation of non-methoxy pentanoic acid-PEG-maleimide (20 kDa)...

SKU: PBSA-01P

Price: \$500.00

Categories: [PEG Proteins](#), [PEG Reagents](#)




MONO-PEGYLATED BSA (MPEG 5 KDA), PBSA-02, 0.1 MG

Prepared by conjugation of mPEG-maleimide (5 kDa) to native BSA followed by...

SKU: PBSA-02

Price: \$500.00

Categories: [PEG Proteins](#), [PEG Reagents](#)




MONO-PEGYLATED BSA (MPEG 20 KDA), PBSA-01, 0.1 MG

Prepared by conjugation of mPEG-maleimide (20 kDa) to the single reactive...

SKU: PBSA-01

Price: \$500.00

Categories: [PEG Reagents](#), [PEG Proteins](#)



MULTI-PEGYLATED BSA, PBSA-00, 1 MG

Prepared by reaction of BSA with a 17-fold molar excess of 20 kDa mPEG...

SKU: PBSA-00

Price: \$400.00

Categories: [PEG Proteins](#), [PEG Reagents](#)



MULTI-PEGYLATED KLH, PKLH-01, 1 MG

Prepared by reaction of 20 kDa
mPEG succinimidyl valerate
(mPEG-SVA) with...

SKU: PKLH-01

Price: \$400.00

Categories: [PEG Proteins](#), [PEG Reagents](#)

INDEX

M

- Mono-PEGylated BSA (mPEG 20 kDa) Biotin Conjugate 1
- Mono-PEGylated BSA (mPEG 20 kDa), PBSA-01, 0.1 mg 1
- Mono-PEGylated BSA (mPEG 5 kDa), PBSA-02, 0.1 mg 1
- Mono-PEGylated BSA (pentanoic acid-PEG, 20 kDa), PBSA-01P, 0.1 mg 1
- Multi-PEGylated BSA, PBSA-00, 1 mg 1
- Multi-PEGylated BSA, PBSA-00-550, 1 mg 1
- Multi-PEGylated KLH, PKLH-01, 1 mg 2