



Phone: 610-431-7707
Fax: 610-431-7818
Email: info@lifediagnostics.com
www.lifediagnostics.com

Life Diagnostics, Inc.
P.O. Box 5205
West Chester PA, 19380
USA

Recombinant Human Cardiac Troponin-I

Catalog No.	1220
Form	Liquid
Concentration	1 - 5 mg/ml (by SPARCL™ immunoassay using pure native human cardiac troponin-I as standard)
Analysis	SDS PAGE: Major 24 kDa band
Purity	>90%
Buffer	8 M urea, 50 mM Tris.HCl, 100 mM NaCl, 15 mM beta-mercaptoethanol, 1 mM EDTA, pH 8.0
Storage	-70°C. Avoid repeated freeze-thaws
Source	Recombinant human cardiac troponin-I expressed in E. coli without tags. Purified by troponin-C affinity chromatography.
Comments	Suitable for use in preparation of calibrators and controls.
Manufactured	In the USA