

Horse C-Reactive Protein (CRP)

Catalog No.	8108
Purity	≥90%
Analysis	SDS PAGE: ≈23 kDa band on Laemmli reducing gels
Buffer	50 mM Tris.Cl (pH 7.5), 200 mM NaCl, 2 mM CaCl ₂ , 0.09% NaN ₃ .
Concentration	≥0.25 mg/ml (assuming Abs 0.1% = 1.809)
Storage	≤-20°C
Source	Horse serum
Manufactured	USA

Horse C-Reactive protein (CRP) is a pentameric serum protein belonging to the pentraxin family that is comprised of five ≈23 kDa subunits. It is a positive acute phase reactant, the levels of which increase in response to inflammation, tissue injury or disease. CRP is purified at Life Diagnostics by multiple affinity chromatography steps.