

Chicken Ex-FABP

Catalog No.	EXFABP-5-0.05
Purity	~90%
Size	50 µg
Analysis	SDS PAGE: 20 kDa band under reducing conditions
Buffer	PBS containing 0.09% NaN ₃ , pH 7.2
Concentration	1-2 mg/ml (Abs 0.1% = 1.153)
Storage	-20°C to -80°C. Avoid repeated freeze-thaws
Source	E. coli

Ex-FABP (extracellular fatty acid binding protein) is a positive acute phase protein in chickens. Serum levels increase in response to inflammation and infection. Ex-FABP is expressed in E. coli with a C-terminal 6-His tag.