

Dog Gastrotropin

Catalog No.	GTP-4-0.05
Purity	~90%
Size	50 µg
Analysis	SDS PAGE: ~15.4 kDa band under reducing conditions
Buffer	PBS, 10% glycerol, pH 7.4
Concentration	1-2 mg/ml (determined by Bradford)
Storage	-20°C to -80°C. Avoid repeated freeze-thaws
Source	E. coli

Gastrotropin is a cytosolic fatty acid binding protein that facilitates bile acid transport and recycling in the intestine. In humans it is a potential cancer biomarker. It is expressed in E. coli with an N-terminal 6-His tag.