



Life Diagnostics, Inc.

Phone: 610-431-7707

Fax: 610-431-7818

Email: info@lifediagnostics.com

www.lifediagnostics.com

Life Diagnostics, Inc.

124 Turner Lane

West Chester PA, 19380

USA

Pig Haptoglobin Monoclonal Antibody

Catalog/Clone	HAPT-9-11D10
Description	Mouse monoclonal Anti-pig haptoglobin Purified from serum-free culture medium Isotype: Mouse IgG1-kappa
Purity	>90% by SDS PAGE
Form	Liquid
Concentration	Lot dependent (E1% at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.10% NaN ₃
Storage	2-8°C
Application	ELISA. Because of the oligomeric nature of haptoglobin the antibody may be used as both capture and detection antibody. The antibody does not work in Western blots.
Immunogen	Pure pig haptoglobin
Manufactured	USA