



Life Diagnostics, Inc.

Phone: 610-431-7707  
Fax: 610-431-7818  
Email: [info@lifediagnostics.com](mailto:info@lifediagnostics.com)  
[www.lifediagnostics.com](http://www.lifediagnostics.com)

Life Diagnostics, Inc.  
P.O. Box 5205  
West Chester PA, 19380  
USA

## Chicken I-FABP

Catalog No.	IFABP-5-0.05
Purity	~90%
Size	50 µg
Analysis	SDS PAGE: ~15.9 kDa band under reducing conditions
Buffer	50 mM Tris-HCl, 150 mM NaCl, 10% glycerol, pH 8.0
Concentration	1-2 mg/ml (determined by Bradford)
Storage	-20°C to -80°C. Avoid repeated freeze-thaws
Source	E. coli

I-FABP (intestinal fatty acid binding protein) is a cytosolic protein involved in fatty acid transport in intestinal epithelial cells. It is a potential biomarker of intestinal lesions caused by coccidiosis. It is expressed in E. coli with an N-terminal 6-His tag.