



Phone: 610-431-7707
Fax: 610-431-7818
Email: info@lifediagnostics.com
www.lifediagnostics.com

Life Diagnostics, Inc.
P.O. Box 5205
West Chester PA, 19380
USA

Multi-PEGylated KLH (5 kDa mPEG)

Catalog No.	PKLH-02
Form	Liquid
Concentration	1-3 mg/ml (Assuming an absorbance of 0.295 at 345 nm for a 1 mg/ml solution). Please note: The concentration above refers only to the KLH polypeptide concentration. It does not include the mass contributed by the 5 kDa mPEG chains
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.09% NaN ₃ , pH 7.2
Storage	4°C Short Term, -20°C to -80°C Long Term
Comments	Prepared by reaction of a 120-fold molar excess of 5 kDa mPEG succinimidyl valerate (5 kDa mPEG-SVA) with purified KLH. Free-PEG is removed by gel-filtration
Manufactured	USA